

AVASPEC NEXOS BENCH ONLY SPECTROMETERS

Operation manual



NEED TECHNICAL SUPPORT?

Scan the QR-code or visit www.avantes.com/support
We are happy to help you!

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, by any means, electronic, mechanical, photocopying, recording, or otherwise, without written permission from Avantes BV.

This manual is sold as part of an order and subject to the condition that it shall not, by way of trade or otherwise, be lent, re-sold, hired out, or otherwise circulated without the prior consent of Avantes BV. in any form of binding or cover other than that in which it is published.

Every effort has been made to make this manual as complete and as accurate as possible, but no warranty or fitness is implied. The information provided is on an "as is" basis. Avantes BV. shall have neither liability nor responsibility to any person or entity concerning any loss or damages arising from the information contained in this manual.

Table of Contents

Introduction	3
1.1 Previous versions and revision codes	3
1.2 AvaSpec- NEXOS™ Bench Only introduction	4
1.3 Technical specifications AvaSpec-NEXOS™ Bench Only.....	4
1.4 Interface cables.....	5
1.4.1 Mating parts	5
1.4.2 Pinout ZIF 8.....	5
1.5 EEPROM.....	6
1.6 Temperature sensor	7
1.7 Dimensions AvaSpec-NEXOS™ Bench Only	7
Appendix A: Detector timing specification	8

Introduction

Thank you for purchasing an Avantes optical bench.

This manual provides OEM users with instructions and technical details on implementing the AvaSpec-NEXOS™ Bench Only into their system.

All information in this manual must be considered confidential and may not be duplicated or given to third parties without written permission from Avantes BV.

1.1 Previous versions and revision codes

OEM Manual – June 2023
OEM Manual – April 2024

First release
Revised product name and Temperature range

1.2 AvaSpec- NEXOS™ Bench Only introduction

The AvaSpec-NEXOS™ Bench Only is a compact optical bench integrated into a small aluminum unibody, as illustrated in Figure 1.

Currently supported detectors are the Hamamatsu S11639-11 (2048CL) and Hamamatsu S13496-11 (4096CL)

The AvaSpec-NEXOS™ Bench Only detector board provides separated video and digital control datalines and contains a 128 kb EEPROM.




Figure 1: AvaSpec-NEXOS™ Bench Only

Specifications regarding detector timing of the control datalines and video signal can be found in Appendix A.

The used interfaces are a ZIF 8 and a coax connector.

The ZIF 8 connector contains the digital control signals and 5VDC, the coaxial connector provides the analog video output signal.

1.3 Technical specifications AvaSpec-NEXOS™ Bench Only

Table 1: Technical specifications AvaSpec-NEXOS™ Bench Only

Digital control	ZIF 8: 8pins, digital signals
Analog ouptut	Coax: 1 Analog video output
Power supply	5VDC, 64 mA The power supply for the AvaSpec-NEXOS™ Bench Only must be able to deliver 5V +/- 5% with a noise level of maximum 20mV. NOTE: The power inlet for the AvaSpec-NEXOS™ Bench Only does not have polarity protection
Temperature range	5 - 55 °C
Dimensions	105 x 80 x 20 mm
Weight	255 g

1.4 Interface cables

The ZIF 8 interface connector is intended to provide digital control signals and power to the sensor. The coax connector is the analog video output of the sensor.

The length of the outgoing flat cable is 130mm (±5mm) with top side contacts. The length of the outgoing coax cable is 195mm (±5mm) with a female coax connector.

1.4.1 Mating parts

Necessary mating parts:

ZIF:

Würth series: FPC Connector and FFC cable 0.5mm SMT ZIF 8 pins horizontal bottom contact.

Würth: 687 108 149 022 - WR-FPC 0.50 mm SMT ZIF Horizontal - Hinge Type

Mini coax connector:

Würth 636 101 111 001 - WR-UMRF PCB Receptacle SMT with 3 Pads

Hirose U.FL-R-SMT(01) male - Light Weight Micro SMT Coaxial Connector, 1.9mm to 2.4mm Mated Height

Mating parts are suggestions, other similar contra connectors can be used.

1.4.2 Pinout ZIF 8

The table and picture below describes the required pinout and orientation of the ZIF 8 cable.

Pinout ZIF 8	Function
	2048CL / 4096CL
1	n.c.
2	Ground
3	Sensor CLK
4	Sensor ST
5	Vdd, 5Volt
6	n.c.
7	I ² C SDA
8	I ² C SCL

Table 2: Interface pinout.




Figure 2: ZIF 8 cable orientation.

1.5 EEPROM

The AvaSpec-NEXOS™ Bench Only is equipped with a non-volatile memory.

As memory the following 128kb(=16kB) serial I²C EEPROM is used:

- 24LC128T-E/MNY from Microchip.

This device uses device address: -0x1010000

For detailed read and write operations please refer to the manufacturer's datasheet.

<http://ww1.microchip.com>

The EEPROM does contain the following information at the given address: see Table 3.

Address	Length	Name	Description
0x00 – 0x01	2 bytes [UINT16]	Structure Version Number	Version number of the data structure in this EEPROM
0x02 – 0x0B	10 bytes	Avabench Serial number	9 characters + 1 terminating zero byte
0x0C – 0x4B	64 bytes	Detector Name	HAM S11639 or HAM S13496
0x4C – 0x4D	2 bytes [UINT16]	Nr total pixels	Total number of data, optical black and invalid pixels
0x4E – 0x4F	2 bytes [UINT16]	Nr obp left	Number of optical black pixels before data pixels, starting at pixel 0. Optical black pixel numbers are: [0 ... Nr_obp_left-1]
0x50 – 0x51	2 bytes [UINT16]	Nr obp right	Number of optical black pixels after data pixels, starting at pixel "Nr_total_pixels -/- Nr_obp_right". Optical black pixel numbers are: [Nr_total_pixels - Nr_obp_right ... Nr_total_pixels-1]
0x52 – 0x53	2 bytes [UINT16]	First effective pixel	First pixel with data
0x54 – 0x55	2 bytes [UINT16]	Last effective pixel	Last pixel with data
0x56 – 0x91	60 bytes 30x [UINT16]	Defective/Hot pixels	Pixels numbers of defective or hot pixels. 30x65535 if no defective or hot pixels
0x92 – 0xA5	20 bytes 5x [FLOAT]	Wavelength calibration	4 th order polynomial for wavelength calibration coefficient (Avantes uses a 3 rd order polynomial (x0 to x3), so x4 = 0.0)

Table 3: EEPROM content.

1.6 Temperature sensor

The new AvaSpec-NEXOS™ Bench Only is equipped with a temperature sensor. This sensor is placed near the detector to obtain the best possible detector temperature. It does not provide a good absolute detector temperature, but can be used to determine temperature changes.

The temperature sensor is accessible via the same I²C bus as the EEPROM, but uses a different device address.

The type of the I²C temperature sensor is:

- LM73CIMK-0/NOPB from Texas Instruments


This temperature sensor uses device address: -0x1001001

For detailed read and write operations please refer to the manufacturer's datasheet.

<https://www.ti.com>


Note: The ALERT pin is not connected, therefore the Alert function is not available.

1.7 Dimensions AvaSpec-NEXOS™ Bench Only




Appendix A: Detector timing specification


source: Hamamatsu datasheet



Detector timing specification Hamamatsu S11639-11 (2048)



Detector timing specification Hamamatsu S13496-11 (4096)



Parameter	Symbol	Min.	Typ.	Max.	Unit
Start pulse width interval *1	tpi(ST)	106/f(CLK)	-	-	s
Start pulse high period *1 *2	thp(ST)	6/f(CLK)	-	-	s
Start pulse low period	tlp(ST)	100/f(CLK)	-	-	s
Start pulse rise and fall times	tr(ST), tf(ST)	0	10	30	ns
Clock pulse duty	-	45	50	55	%
Clock pulse rise and fall times	tr(CLK), tf(CLK)	0	10	30	ns

*1: Dark output increases if the start pulse period or start pulse high period is lengthened.

*2: The integration time equals the high period of ST plus 48 CLK cycles.

The shift register starts operation at the rising edge of CLK immediately after ST goes low.

The integration time can be changed by changing the ratio of the high and low periods of ST.