

Avantes Announces AvaSoft Release 8.9.2.0

Important Software Update

Date: February 27, 2019

Appeldorn, The Netherlands— February 26, 2019 Avantes is excited to announce the release of the AvaSoft 8.9.2.0 software update as well as updates to our AvaSpec DLL and Linux libraries.

“Advanced users will appreciate the new features and novice users will find saving and exporting data more intuitive.” said Ryan Flaherty, Regional Sales Manager with Avantes

Positive Customer Impact

In addition to advances to Ethernet connectivity and stability, this release also provides enhancements to the saving features in experiment files as well as for the periodic saving of merged spectral files. Users will also appreciate the new 3D mode for making custom reflectance measurements using a calibrated standard.

Update Availability

Update is an update driven by customer feedback and is part of Avantes’ commitment to support the scientists and industries that rely on Avantes spectrometers. Update is available for immediately download at <https://www.avantes.com/support/software>. Website registration is required for download.

New DLL and Linux Library

To support a variety of new detectors, for customers that write their own software, Avantes has released the new AvaSpec-DLL version v9.8.1.0. The DLL also includes a Linux library and includes sample programs in Labview, Matlab, C#, C++, Visual Studio, Python QT and Delphi. The DLL is available for 32 and 64 bit operating systems. For more information about obtaining the latest version of the DLL, please contact Avantes at infoUSA@avantes.com or info@avantes.com.

Avantes Enlightening Spectroscopy

Founded in 1994, Avantes the leading innovator in the development of fiber optic spectroscopy instrumentation. Avantes offers a wide range of state of the art, innovative measuring equipment that helps mankind live longer, healthier, and saving the planet for generations to come.

For more information, press only:

Caroline Bach
+31-(0)-313-670170
c.bach@avantes.com

For more information on Update:

<https://www.avantes.com/support/software>

