

Direct-attach Accessories

CUV-DA



Avantes offers a wide range of high quality direct-attach accessories, such as cuvette holders, filter holders and attenuators for the AvaLight series of light sources. In the table below an overview of direct-attach accessories for each light source can be found. For more information, please take a look at the following pages.

The last column is a combination of the attenuator, the cuvette holder and the filter holder in one easy to attach device.

CUV-ATT-DA



Ordering Information

Coupling to	Mounting	Round 1/2" Filter Holder	Cuvette Holder	Attenuator	Cuvette/ Filter Holder/ Attenuator
AvaLight-HAL-Mini	UNS thread	FH-DA-HAL-Mini	CUV-DA-HAL-Mini	ATT-DA-HAL-Mini	CUV-ATT-DA-HAL-Mini
AvaLight-DHS	Front plate	n/a	CUV-DA-DHS	n/a	n/a
AvaLight-DHc					
AvaLight-XE	SMA thread	FH-DA	CUV-DA	ATT-DA	CUV-ATT-DA
AvaLight-LED					
Fiber-optic	fiber	FH-INLINE	CUV-UV/VIS	ATT-INL-EXT	n/a