

Fiber-optic Connectors

SMA Connectors

We supply all of our standard fiber-optic cables, bundles and probes with SMA-905 connectors that easily fit into our complete range of spectrometers, light sources and accessories.

The SMA-905 connectors are screw-fitted and can be rotated over 360 degrees. The typical insertion loss for the connectors is 0.5 dB. The maximum filling diameter for bundles is 2.46 mm

FC/PC Connectors

Optional FC/PC-connectors can be mounted to our fiber-optic products. The multi-mode FC/PC connectors have an extremely low insertion loss of < 0.2 dB. The FC/PC connector cannot rotate, always mounts into the same fixed position and therefore has a high reproducibility.

Ordering Information

- FC/PC • FC/PC connector instead of standard SMA

